

Personal Sharps Container Illustrative User Guide




SAFE NEEDLE DISPOSAL INSTRUCTIONS:

1 Keep new, unused needles in the larger compartment on the left (10). They will be separated from your used needles by a flap. Do not recap used needles. Do not snap, break or bend any part of the needle. **2** Carefully push used needles point-first past the flap on the right. Talk to your distributor about where to dispose of a fully used sharps container.



NEW



- 1.** STORE NEW NEEDLES HERE. TAKE NEW NEEDLE FOR SINGLE USE

↑
OPEN HERE



USED



- 2.** PUSH USED NEEDLE THROUGH FLAP (THIS SIDE). ASK WHERE TO DISPOSE OF WHEN CONTAINER IS FULL (10 NEEDLES)