

HIV Point of Care Test (INSTI™) Quality Control Log

Testing site: _____

Date of Quality Control dd/mm/yyyy	Tester Initials	Quality Control Material Negative, Positive HIV-1	Quality Control Lot Number (On Vial)	POC Test Lot Number (On Back of Kit Package)	POC Test Result Reactive Non-reactive Invalid	POC Test performs as expected (Yes, No)	Notes / Comments Reason for testing (from numbers below)	Site Lead Review Initial

1. Receive new shipment

3. New tester, before they start testing

5. Change to new lot number

7. Temperature out of range (below 2°C or above 30°C)

- 2. Invalid client result X 2
- 4. At regular intervals monthly
- 6. Trouble-shooting