



CUN-18-001

Customer Update Notice

Date: 13-Aug-2018

Re: INSTI® HIV-1/HIV-2 Antibody Test Kit.

Dear Customer,

bioLytical has received approval from regulatory authorities for extending temperature storage from 15-30°C (59 – 86°F) to 2-30°C (35.6 – 86°F) for INSTI® HIV-1/HIV-2 Antibody Test Kit.

Based on this customer update notice, you can store the kits at 2-30°C (35.6 – 86°F).

Please note that for a limited period, you may still receive product with an older version of the membrane unit pouch labeling, kit pouch labeling, or Instruction for Use (IFU) with 15-30°C storage temperature until the existing inventory of the labels and IFUs are used.

Please contact bioLytical Laboratories Inc. customer support if you have any questions or concerns at +1-866-674-6784, or by email at customercare@biolytical.com.

Sincerely,

A handwritten signature in black ink, appearing to read 'Hans Croukamp', written over a circular stamp or mark.

Hans Croukamp

Senior Director, Operations